

# CRMP5 antibody

Catalog Number: 109574

## Product name

CRMP5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CRMP5 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IHC detected in human prostate cancer tissue. Observed molecular weight by Western-blot: 61kd

## Preparation

This antibody was obtained by immunization of CRMP5 recombinant protein (Accession Number: NM\_020134). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109574(DPYSL5 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:109574(DPYSL5 antibody) at dilution of 1:100 (under 25x lens)