

# CRYBA2 antibody



Catalog Number: 109583

## Product name

CRYBA2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CRYBA2 Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue, rat eye tissue. Observed molecular weight by Western-blot: 22 kDa

## Preparation

This antibody was obtained by immunization of CRYBA2 recombinant protein (Accession Number: NM\_057093). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:109583(CRYBA2 Antibody) at dilution of 1:2000