

CRYBB3 antibody



Catalog Number: 109585

Product name

CRYBB3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CRYBB3 Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue. Positive IF detected in HepG2 cells, Hela cells. Observed molecular weight by Western-blot: 24 kDa

Preparation

This antibody was obtained by immunization of CRYBB3 recombinant protein (Accession Number: XM_047441147). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:109585(CRYBB3 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using CRYBB3 antibody Catalog No:109585 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).