

## Product name

CSN3 antibody

## Specificity

Human; other species not tested.

## Antibody description

CSN3 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue.

## Preparation

This antibody was obtained by immunization of CSN3 recombinant protein (Accession Number: XM\_011531620). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:10-1:100

## Validations



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:109599(CSN3 Antibody) at dilution of 1:25 (under 10x lens)



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:109599(CSN3 Antibody) at dilution of 1:25 (under 40x lens)