

# CSNK1A1 antibody



Catalog Number: 109600

## Product name

CSNK1A1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CSNK1A1 Rabbit Polyclonal antibody. Positive IF detected in A549 cells. Positive WB detected in MCF7 cells, HeLa cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 42-45 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001892).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:109600(CSNK1A1 antibody) at dilution of



Immunofluorescent analysis of A549 cells using Catalog No:109600(CSNK1A1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)