

c-SRC antibody



Catalog Number: 109607

Product name

Polyclonal

c-SRC antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

c-SRC Rabbit Polyclonal antibody. Positive WB detected in A431 cells, CHO cells, COLO 320 cells, Jurkat cells, mouse brain tissue. Observed molecular weight by Western-blot: 60 kDa

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of c-SRC recombinant protein (Accession Number: XM_047440405). Purification method: Antigen affinity purified.

Formulation


PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:109607(c-SRC Antibody) at dilution of 1:1000