

# CST8 antibody



Catalog Number: 109614

## Product name

CST8 antibody

## Specificity

Human; other species not tested.

## Antibody description

CST8 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, DU 145 cells, HEK-293 cells. Observed molecular weight by Western-blot: 17-19 kDa, 9 kDa

## Preparation

This antibody was obtained by immunization of CST8 recombinant protein (Accession Number: XM\_047439819). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109614(CST8 antibody) at dilution of 1:1000