

# CTAGE1 antibody



Catalog Number: 109621

## Product name

CTAGE1 antibody

## Specificity

Human; other species not tested.

## Antibody description

CTAGE1 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 86-120 kDa

## Preparation

This antibody was obtained by immunization of CTAGE1 recombinant protein (Accession Number: BC031065). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109621(CTAGE1 antibody) at dilution of 1:800



IP Result of anti-CTAGE1 (IP:Catalog No:109621, 4ug; Detection:Catalog No:109621 1:500) with HeLa cells lysate 2400ug.