

CTAGE6 antibody



Catalog Number: 109623

Product name

CTAGE6 antibody

Specificity

Human; other species not tested.

Antibody description

CTAGE6 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, MCF-7 cells. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of CTAGE6 recombinant protein (Accession Number: BC043153). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109623(CTAGE6 Antibody) at dilution of 1:600