

CTAGE6 antibody



Catalog Number: 109623

Product name

Polyclonal

CTAGE6 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

CTAGE6 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, MCF-7 cells. Observed molecular weight by Western-blot: 100 kDa

Dilutions

Recommended Dilution:


Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of CTAGE6 recombinant protein (Accession Number: BC043153). Purification method: Antigen affinity purified.

Validations

Formulation



PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109623(CTAGE6 Antibody) at dilution of 1:600

Store at -20°C. DO NOT ALIQUOT

Clonality