

# CXADR antibody

Catalog Number: 109635

## Product name

CXADR antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CXADR Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in BxPC-3 cells, mouse heart tissue, mouse small intestine tissue. Observed molecular weight by Western-blot: 55 kDa

## Preparation

This antibody was obtained by immunization of CXADR recombinant protein (Accession Number: NM\_001338). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

## Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109635(CXADR antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:109635(CXADR antibody) at dilution of 1:50 (under 10x lens)