

CXCL13/BCA1 antibody



Catalog Number: 109638

Product name

CXCL13/BCA1 antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

CXCL13/BCA1 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human breast cancer tissue, human liver tissue, human lymphoma tissue. Positive WB detected in human placenta tissue, human serum tissue, Jurkat cells. Observed molecular weight by Western-blot: 9 kDa

Preparation

This antibody was obtained by immunization of CXCL13/BCA1 recombinant protein (Accession Number: NM_001371558). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:109638(CXCL13 antibody) at dilution of 1:200



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109638(CXCL13 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded

CXCL13/BCA1 antibody



Catalog Number: 109638

human tonsillitis tissue slide using Catalog

No:109638(CXCL13 Antibody) at dilution of 1:50
(under 40x lens)