

CYB5R4 antibody



Catalog Number: 109658

Product name

Polyclonal

CYB5R4 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

CYB5R4 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, HeLa cells. Observed molecular weight by Western-blot: 55.8kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of CYB5R4 recombinant protein (Accession Number: NM_016230). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109658(CYB5R4 antibody) at dilution of 1:500

