

# CYB5R4 antibody



Catalog Number: 109658

## Product name

CYB5R4 antibody

## Immunogen

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

CYB5R4 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, HeLa cells. Observed molecular weight by Western-blot: 55.8kd

## Preparation

This antibody was obtained by immunization of CYB5R4 recombinant protein (Accession Number: NM\_016230). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109658(CYB5R4 antibody) at dilution of 1:500