

Cyclin T1 antibody



Catalog Number: 109670

Product name

Cyclin T1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Cyclin T1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, A431 cells, HeLa cells, Jurkat cells, Y79 cells. Observed molecular weight by Western-blot: 87 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001240). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:109670(CCNT1 antibody) at dilution of 1:300