

# Cyclophilin E antibody



Catalog Number: 109675

## Product name

Polyclonal

Cyclophilin E antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

Cyclophilin E Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 33 kDa

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of Cyclophilin E recombinant protein (Accession Number: NM\_006112). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109675(PPIE Antibody) at dilution of 1:1500