

CYFIP2 antibody

Catalog Number: 109676

Product name

CYFIP2 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

CYFIP2 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, mouse brain tissue. Positive IHC detected in human gliomas tissue, human brain tissue. Observed molecular weight by Western-blot: 130 kDa

Preparation

This antibody was obtained by immunization of CYFIP2 recombinant protein (Accession Number: XM_011534516). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109676(CYFIP2 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:109676(CYFIP2 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human gliomas using Catalog No:109676(CYFIP2 antibody) at dilution of 1:50 (under 10x lens)

CYFIP2 antibody

Catalog Number: 109676

antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-CYFIP2 (IP:Catalog No:109676, 4ug; Detection:Catalog No:109676 1:700) with mouse brain tissue lysate 4000ug.