

CYP20A1 antibody



Catalog Number: 109686

Product name

CYP20A1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CYP20A1 Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 45kd,60kd

Preparation

This antibody was obtained by immunization of CYP20A1 recombinant protein (Accession Number: NM_177538). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109686(CYP20A1 antibody) at dilution of 1:2000



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:109686(CYP20A1 antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:109686(CYP20A1 antibody) at dilution of 1:200 (under 40x lens)