

CYP24A1 antibody



Catalog Number: 109687

Product name

CYP24A1 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CYP24A1 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, A549 cells. Positive IHC detected in human kidney tissue, human prostate cancer tissue. Observed molecular weight by Western-blot: 54-59 kDa

Preparation

This antibody was obtained by immunization of CYP24A1 recombinant protein (Accession Number: XM_047439936). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109687(CYP24A1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human kidney using Catalog No:109687(CYP24A1 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human kidney using Catalog No:109687(CYP24A1 antibody) at dilution of 1:50 (under 10x lens)

CYP24A1 antibody



Catalog Number: 109687

antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-CYP24A1 (IP:Catalog No:109687, 4ug; Detection:Catalog No:109687 1:300) with HeLa cells lysate 2800ug.