

# CYP2F1 antibody



Catalog Number: 109695

## Product name

CYP2F1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

CYP2F1 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 50-60kd

## Preparation

This antibody was obtained by immunization of CYP2F1 recombinant protein (Accession Number: XM\_047438282). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:109695(CYP2F1 antibody) at dilution of 1:500