

CYP2S1 antibody



Catalog Number: 109697

Product name

CYP2S1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CYP2S1 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, mouse kidney tissue. Positive IF detected in Hela cells. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of CYP2S1 recombinant protein (Accession Number: NM_030622). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:109697(CYP2S1 antibody) at dilution of 1:500



Immunofluorescent analysis of Hela cells, using CYP2S1 antibody Catalog No:109697 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).