

CYP4F11 antibody



Catalog Number: 109705

Product name

CYP4F11 antibody

IHC: 1:20-1:200

Specificity

human; other species not tested.

Validations

Antibody description

CYP4F11 Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer. Positive WB detected in human liver tissue, human kidney tissue. Observed molecular weight by Western-blot: 57kd




Preparation

This antibody was obtained by immunization of CYP4F11 recombinant protein (Accession Number: NM_021187). Purification method: Antigen affinity purified.

human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109705(CYP4F11 antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:109705(CYP4F11 antibody) at dilution of 1:100 (under 10x lens)

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunohistochemical of paraffin-embedded human liver cancer using Catalog

CYP4F11 antibody



Catalog Number: 109705

No:109705(CYP4F11 antibody) at dilution of 1:100 (under 40x lens)