

# CYP4F12 antibody



Catalog Number: 109707

## Product name

CYP4F12 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CYP4F12 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 60kd

## Preparation

This antibody was obtained by immunization of CYP4F12 recombinant protein (Accession Number: NM\_023944). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:109707(CYP4F12 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109707(CYP4F12 antibody) at dilution of 1:100 (under 40x lens)