

CYP4Z1-Specific antibody



Catalog Number: 109710

Product name

CYP4Z1-Specific antibody

Specificity

Human; other species not tested.

Antibody description

CYP4Z1-Specific Rabbit Polyclonal antibody.
Positive WB detected in MDA-MB-453s cells, HeLa cells. Observed molecular weight by Western-blot: 59kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_178134).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:109710(CYP4Z1-Specific antibody) at dilution of 1:300