

Cystatin A antibody



Catalog Number: 109714

Product name

Cystatin A antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Cystatin A Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse liver tissue, A549 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 11 kDa

Preparation

This antibody was obtained by immunization of Cystatin A recombinant protein (Accession Number: NM_005213). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109714(CSTA antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using CSTA antibody Catalog No:109714 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).