

# CTSO antibody



Catalog Number: 109735

## Product name

CTSO antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CTSO Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 36 kDa

## Preparation

This antibody was obtained by immunization of CTSO recombinant protein (Accession Number: NM\_001334). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109735(CTSO Antibody) at dilution of 1:700