

CUL4A antibody

Catalog Number: 109744

Product name

CUL4A antibody

Specificity

Human; other species not tested.

Antibody description

CUL4A Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue.

Preparation

This antibody was obtained by immunization of CUL4A recombinant protein (Accession Number: NM_001008895). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:109744(CUL4A antibody) at dilution of 1:100 (under 10x lens)