

# DCTN4 antibody



Catalog Number: 109759

## Product name

DCTN4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DCTN4 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in human brain tissue, HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 62 kDa

## Preparation

This antibody was obtained by immunization of DCTN4 recombinant protein (Accession Number: NM\_016221). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109759(DCTN4 antibody) at dilution of 1:800



IP Result of anti-DCTN4 (IP:Catalog No:109759, 4ug; Detection:Catalog No:109759 1:700) with Jurkat cells lysate 2800ug.