

DCTPP1 antibody

Catalog Number: 109762

Product name

DCTPP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DCTPP1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse heart tissue, Transfected HEK-293 cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 19 kDa; 40 kDa; 80 kDa

Preparation

This antibody was obtained by immunization of DCTPP1 recombinant protein (Accession Number: NM_024096). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109762(DCTPP1 Antibody) at dilution of 1:300



Immunofluorescent analysis of HeLa cells using Catalog No:109762(DCTPP1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)