

# DCX antibody

Catalog Number: 109766

## Product name

DCX antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DCX Rabbit Polyclonal antibody. Positive IHC detected in human ovary tissue, human placenta tissue, human spleen tissue. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 40kd

## Preparation

This antibody was obtained by immunization of DCX recombinant protein (Accession Number: NM\_178151). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109766(DCX antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human ovary using Catalog No:109766(DCX antibody) at dilution of 1:50 (under 10x lens)