

DDI1 antibody

Catalog Number: 109773

Product name

DDI1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DDI1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, U-937 cells. Observed molecular weight by Western-blot: 47kDa, 56 kDa

Preparation

This antibody was obtained by immunization of DDI1 recombinant protein (Accession Number: NM_001001711). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109773(DDI1 antibody) at dilution of 1:1000