

# DDO antibody

Catalog Number: 109774

## Product name

DDO antibody

## Immunogen

[Human DDO Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DDO Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IP detected in SKOV-3 cells. Observed molecular weight by Western-blot: 40 kDa

## Preparation

This antibody was obtained by immunization of DDO recombinant protein (Accession Number: NM\_003649). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109774(DDO Antibody) at dilution of 1:700



IP Result of anti-DDO (IP:Catalog No:109774, 4ug; Detection:Catalog No:109774 1:500) with SKOV-3 cells lysate 2800ug.