

CYTB antibody



Catalog Number: 109782

Product name

CYTB antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CYTB Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in MCF-7 cells. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: YP_003024038). Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:109782(CYTB Antibody) at dilution of 1:300



Immunofluorescent analysis of HeLa cells using Catalog No:109782(CYTB Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)