

CYTSB antibody



Catalog Number: 109816

Product name

CYTSB antibody

Immunogen

[Human CYTSB Recombinant protein \(GST tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

CYTSB Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HeLa cells. Observed molecular weight by Western-blot: 119kd

Preparation

This antibody was obtained by immunization of CYTSB recombinant protein (Accession Number: NM_001243438). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:109816(CYTSB antibody) at dilution of 1:1000