

DDX24 antibody



Catalog Number: 109823

Product name

DDX24 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DDX24 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Y79 cells. Observed molecular weight by Western-blot: 120 kDa

Preparation

This antibody was obtained by immunization of DDX24 recombinant protein (Accession Number: NM_020414). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109823(DDX24 antibody) at dilution of 1:300