

DDX25 antibody

Catalog Number: 109824

Product name

DDX25 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

DDX25 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human testis tissue, Jurkat cells. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 50-55 kDa

Preparation

This antibody was obtained by immunization of DDX25 recombinant protein (Accession Number: NM_001330438). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109824(DDX25 antibody) at dilution of 1:500



IP Result of anti-DDX25 (IP:Catalog No:109824, 4ug; Detection:Catalog No:109824 1:500) with mouse testis tissue lysate 4000ug.