

# DDX27 antibody



Catalog Number: 109825

---

## Product name

DDX27 antibody

## Immunogen

[Human DDX27 Recombinant protein \(His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

DDX27 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse liver tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 90kd

## Preparation

This antibody was obtained by immunization of DDX27 recombinant protein (Accession Number: NM\_017895). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109825(DDX27 antibody) at dilution of 1:500