

DEDD2 antibody



Catalog Number: 109847

Product name

DEDD2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DEDD2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 36kd

Preparation

This antibody was obtained by immunization of DEDD2 recombinant protein (Accession Number: NM_133328). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109847(DEDD2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human brain using Catalog No:109847(DEDD2 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:109847(DEDD2 antibody) at dilution of 1:100 (under 40x lens)