

DEFA1 antibody

Catalog Number: 109849

Product name

DEFA1 antibody

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

Human; other species not tested.

Validations

Antibody description

DEFA1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Positive IF detected in HepG2 cells. Positive IHC detected in human small intestine tissue, human liver cancer tissue. Observed molecular weight by Western-blot: 6 kDa




Preparation

This antibody was obtained by immunization of DEFA1 recombinant protein (Accession Number: NM_004084). Purification method: Antigen affinity purified.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109849(DEFA1 Antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:109849(DEFA1 Antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF



Dilutions

Recommended Dilution:


WB: 1:200-1:1000

Immunohistochemistry of paraffin-embedded

DEFA1 antibody

Catalog Number: 109849

human small intestine tissue slide using Catalog No:109849(DEFA1 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:109849(DEFA1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)