

DAB2IP antibody

Catalog Number: 109852

Product name

DAB2IP antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

DAB2IP Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 118 kDa

Preparation

This antibody was obtained by immunization of DAB2IP recombinant protein (Accession Number: BC146762). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109852(DAB2IP Antibody) at dilution of 1:1000



IP Result of anti-DAB2IP (IP:Catalog No:109852, 4ug; Detection:Catalog No:109852 1:1000) with mouse brain tissue lysate 4000ug.