

# DAND5 antibody

Catalog Number: 109857

## Product name

DAND5 antibody

## Specificity

Human; other species not tested.

## Antibody description

DAND5 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Positive IHC detected in human testis tissue, human heart tissue. Observed molecular weight by Western-blot: 20 kDa

## Preparation

This antibody was obtained by immunization of DAND5 recombinant protein (Accession Number: NM\_152654). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109857(DAND5 antibody) at dilution of 1:400



Immunohistochemistry of paraffin-embedded human testis slide using Catalog No:109857(DAND5 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human testis slide using Catalog No:109857(DAND5 Antibody) at dilution of 1:50