

DAP13 antibody



Catalog Number: 109859

Product name

Polyclonal

DAP13 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

DAP13 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 20kd

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of DAP13 recombinant protein (Accession Number: NM_018838). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:109859(NDUFA12 antibody) at dilution of 1:800