

# DAZ4 antibody



Catalog Number: 109872

## Product name

DAZ4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DAZ4 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human testis tissue, rat testis tissue. Observed molecular weight by Western-blot: 55-60 kDa

## Preparation

This antibody was obtained by immunization of DAZ4 recombinant protein (Accession Number: NM\_020420). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109872(DAZ4 Antibody) at dilution of 1:600