

DAZAP1 antibody



Catalog Number: 109873

Product name

DAZAP1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

DAZAP1 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer.

Preparation

This antibody was obtained by immunization of DAZAP1 recombinant protein (Accession Number: NM_018959). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:109873(DAZAP1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:109873(DAZAP1 antibody) at dilution of 1:100 (under 40x lens)