

# DBF4 antibody



Catalog Number: 109876

## Product name

DBF4 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

DBF4 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 70-80kd

## Preparation

This antibody was obtained by immunization of DBF4 recombinant protein (Accession Number: NM\_006716). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109876(DBF4 antibody) at dilution of 1:300