

# DBP antibody



Catalog Number: 109879

## Product name

Polyclonal

DBP antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

DBP Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 35 kDa

## Dilutions

Recommended Dilution:


## Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of DBP recombinant protein (Accession Number: NM\_001352). Purification method: Antigen affinity purified.

## Validations

## Formulation



PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:109879(DBP Antibody) at dilution of 1:300

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality