

DGCR2 antibody



Catalog Number: 109904

Product name

DGCR2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DGCR2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 61-65 kDa

Preparation

This antibody was obtained by immunization of DGCR2 recombinant protein (Accession Number: NM_005137). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109904(DGCR2 antibody) at dilution of 1:300