

dGK antibody

Catalog Number: 109907

Product name

dGK antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

dGK Rabbit Polyclonal antibody. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of dGK recombinant protein (Accession Number: NM_080916). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-DGUOK (IP:Catalog No:109907, 3ug; Detection:Catalog No:109907 1:300) with mouse kidney tissue lysate 4000ug.