

# DHRS1 antibody

Catalog Number: 109921

## Product name

DHRS1 antibody

## Immunogen

[Human DHRS1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DHRS1 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, human liver tissue, mouse testis tissue. Positive IHC detected in human kidney tissue, human heart tissue, human liver tissue, human ovary tissue, human skin tissue, human testis tissue. Observed molecular weight by Western-blot: 35-38 kDa

## Preparation

This antibody was obtained by immunization of DHRS1 recombinant protein (Accession Number: NM\_001136050). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:109921(DHRS1 antibody) at dilution of 1:800




Immunohistochemical of paraffin-embedded human kidney using Catalog No:109921(DHRS1 antibody) at dilution of 1:50 (under 10x lens)

# DHRS1 antibody

Catalog Number: 109921

---



Immunohistochemical of paraffin-embedded human kidney using Catalog No:109921(DHRS1 antibody) at dilution of 1:50 (under 40x lens)