

DICER1 antibody

Catalog Number: 109939

Product name

DICER1 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

DICER1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human lung cancer tissue. Positive IF detected in HepG2 cells. Positive IP detected in Jurkat cells. Positive WB detected in K-562 cells, Jurkat cells. Observed molecular weight by Western-blot: 220-250kd, 90kd

Validations



Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_030621). Purification method: Antigen affinity purified.

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:109939(DICER1 antibody) at dilution of 1:200

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF, IP


Dilutions

Recommended Dilution:

WB: 1:200-1:1000




IP Result of anti-DICER1 (IP:Catalog No:109939, 5ug; Detection:Catalog No:109939 1:300) with Jurkat cells lysate 2000ug.



DICER1 antibody

Catalog Number: 109939

Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:109939(DICER1 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:109939(DICER1 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG